Betriebsanleitung Temperatur-Kalibrierset
Operating instruction temperature calibrating set
Mode d'emploi jeu de calibrage de la température

KERN DBS-A01

Version 1.1

02/2013

Deutsch

English

Français

DBS_A01-BA-def-1311

KERN DBS-A01

Version 1.1 02/2013

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2 Appliance overview

2.1 Overview of display

Display	Description
MAX	Means the highest measured value of a measurement
MIN	Means the lowest measured value of a measurement
HOLD	Means a recorded display value
HI OFF	Alarm for the upper temperature limit value switched off
LO OFF	Alarm for the lower temperature limit value switched off
HI XXX	Upper temperature limit value
Lo XXX	Lower temperature limit value

2.2 Keyboard and function overview

	Description
ON OFF	Switch appliance on/off
	 During display of a MAX/MIN or a HOLD value the unit cannot be switched off.
UP	Decrease value
HOLD	Press the button to record the current displayed value, the indicator „HOLD“ flashes. It remains displayed until it is deleted by pressing button again.

DOWN	Increase value
MAX/MIN	Display Max. / Min. temperature Press button, the highest measured value will be displayed, the indicator „MAX“ flashes. Press button again, the lowest measured value will be displayed, the indicator „MIN“ flashes. Press button anew, the appliance returns to measuring mode. To delete the MAX/MIN values, switch appliance off and on using the ON/OFF feature.
SET	Setting alarm for temperature limit values Press [SET], „HI XXX“ flashes. St the upper limit value using [UP] and [DOWN]. Every time the button is pressed, the value is increased/reduced. If the button is permanently, the value is automatically increased/reduced. Use [ON/OFF] to switch off the limit value, „HI OFF“ is displayed. Confirm upper limit value using [SET], „Lo XXX“ flashes. Adjust the lower limit value in the some way as described at the upper limit value. If the temperature is beyond the limit values, the alarm will be triggered. The alarm can be switched off with any button. The indicator [HI] or [LO], however, remains as long as the temperature is beyond the limit values. Commutation °C / °F Keep set [SET] pressed, the display changes from °C ⇔ °F.

3 Change battery

If the battery tension is not more sufficient, the display darkens.

- ☞ Turn off unit.
- ☞ Open the lid of the battery compartment loosening the 4 screws.
- ☞ Replace battery.

Ensure correct polarisation.


- ☞ Close again the lid.

4 Adjust temperature value of the moisture analyzer KERN DBS


We recommend sometimes to check the temperature value of the appliance using the optional temperature calibrating set DBS-A01. The temperature is measured at two points (100°C and 180°C) and can be adjusted at these two points. The equipment should cool down first to room temperature after the last heating period.

Preparation:


- ⇒ Fix the temperature sensor on the temperature calibration set in accordance with the illustration.


- ⇒ Remove sample dish and dish holder from the moisture analyser.

- ⇒ Install the temperature-calibration set acc. to fig.


- ⇒ Close heated cover.
If applicable, the error message "ol" will be displayed, but can be ignored in this case.
- ⇒ On the temperature calibration set switch-on the digital thermometer by the **ON/OFF** button.

Call up menu:

- ⇒ Press the **Menu** button to access the menu and the first menu item "**PRoGRM**" will be displayed.
- ⇒ Use navigation buttons **↓** **↑** to select menu item "**CAL**".
- ⇒ Acknowledge using **ENTER** key, the current setting is displayed. If necessary, use navigation buttons **↓** **↑** to select "**TEMP**".
- ⇒ Press **ENTER** to confirm and the first temperature point to be tested will be displayed.
- ⇒ **Press START** and the first heating-up period will start.

- Do not leave the heated cover open for more than 1 minute during the heating-up period, otherwise the error message "**ERR.100**" will appear. In this case press **ESC** and restart menu item "**TEMP**".
- Attention: some parts e.g. ventilation grilles and inspection windows may become very hot during operation.
- The adjustment must be made within 15 minutes when the temperature display flashes, otherwise the temperature adjustment will be aborted (AboRT will be displayed).

- ⇒ Temperature calibration for the first point takes 15 min. Compare the displayed temperature value on the thermometer with that of the moisture analyser whilst the temperature display is flashing. If the two values do not match, correct them using the navigation buttons **↓** **↑** and confirm with the **ENTER** button. Numeric input see operating instructions KERN DBS, chap. 7.1.1.

⇒ The second temperature point to be tested will be displayed.

⇒ Confirm with the **ENTER** button, the second heat-up phase is started.

⇒ Temperature calibration for the second point takes 15 min. Compare the displayed temperature value on the thermometer with that of the moisture analyser whilst the temperature display is flashing. If the two values do not match, correct them using the navigation buttons **↓** **↑** and confirm with the **ENTER** button. Numeric input see operating instructions KERN DBS, chap. 7.1.1.

⇒ The adjustment is completed when “END” is displayed. The appliance returns to the menu automatically.
The adjustment protocol output is started upon connection to an optional printer and activated GLP function. (see operating instructions KERN DBS chap. 6.3).

⇒ Press **ESC** to return to Moisture analysis mode.

- The adjustment protocol output is started upon connection to an optional printer and activated GLP function. (see operating instructions KERN DBS chap. 6.3).

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